



# Tucson Python Meetup

How I Use Image Analysis & Processing at Work


Franklin Ventura

March 8, 2016



# Franklin Ventura

- BS Biomedical Engineering UofA
- MS Software/Systems Engr ASU
- Ventana Medical Systems 5 yrs
  - Biomedical/Software Engineer
  - Technology and Applied Research
  - Digital Pathology
  - Cancer Diagnostics
- Applied Imaging Systems Team
  - Karl Garsha, Dustin Harshman, Taras Golota



# Research Development







◎ Flexible System

◎ Fluorescence

○ WS 10-80GB

○ Study ~3TB

“




## Core Libraries

- ◎ VTK
  - Image rendering and manipulation
- ◎ ITK
  - Image registration
- ◎ Matplotlib
  - Plotting
- ◎ Scikit-learn
  - Machine learning
- ◎ Scikit-image
  - Image processing
- ◎ Scipy
  - Scientific functions
- ◎ PIL
  - Simple image display and manipulation
- ◎ Numpy
  - Math functions, array manipulation
- ◎ Omero
  - Image database integration
- ◎ LaTeX/ReportLab
  - Report/Documentation Generation
- ◎ Micromanager
  - Connect to camera, microscope, etc.
- ◎ PywinAuto
  - Automation
- ◎ cx\_Freeze and Innosetup
  - Compiler and installation
- ◎ wxPython
  - GUI
- ◎ mlabwrap
  - Matlab incorporation


# Visualization Toolkit (VTK) - Kitware

*3D computer graphics, modeling,  
image processing, volume  
rendering, scientific visualization,  
and information visualization*





# Omero - Image Database



# Micro-Manager - Hardware Control


## Motorized Microscopes

- **iMic** - From Till Photonics
- **LeicaDMI** - Leica DMI and DM microscopes
- **LeicaDMR** - Leica DMR microscopes (and C
- **MoticMicroscope** - Motic BA600 microscop
- **NikonTE2000** - Nikon TE2000 motorized m
- **NikonTI** - Nikon Ti microscope
- **NikonAZ100** - Nikon AZ100 Zoom microsc
- **Olympus** - Olympus IX81 (Inverted)/BX61
- **ZeissCAN** - Zeiss CAN-bus (not CAN29) co
- **ZeissCAN29** - Zeiss CAN29-bus compatibl

## Cameras


- **ABSCamera** - Cameras from ABS Gesellsc Jena
- **Andor** - Andor cameras
- **AndorSDK3** - New model Andor cameras
- **Apogee** - Apogee cameras
- **AxioCam** - Zeiss AxioCam medium resolut
- **BaumerOptron** - Leica DFC camera (thru
- **dc1394** - Firewire cameras that adhere to i
- **DECamera** - Direct Electron Cameras
- **DemoCamera** - Virtual, emulated equipme
- **DemoStreamCamera** - Virtual, emulated
- **DSLRRemoteCamera** - Canon DSLR came
- **GigECamera** - GigE Vision- and GenICam-
- **GPhoto** - Canon and Nikon DSLR cameras
- **Hamamatsu** - Hamamatsu cameras (throu

```
import MMCorePy
mmc = MMCorePy.CMMCore()
mmc.loadDevice("cam", "PVCAM",
"Camera-1")
mmc.initializeDevice("cam")
mmc.setProperty("cam", "Exposure", 100)
mmc.setProperty("cam", "Binning", 2)
mmc.snapImage()
img = mmc.getImage()
```



## ReportLab


“Create documents in Adobe's Portable Document Format (PDF) using the Python programming language






## cx\_Freeze and Inno Setup

- ◎ cx\_Freeze
  - Cross platform freeze Python scripts into executables
- ◎ Inno Setup
  - Installer for Windows programs






# Scipy and Numpy






## [EX1: Tissue Colorization \(Click for example\)](#)

Hematoxylin and Eosin (H&E) staining is most commonly used stains in histology. The stain enables pathologists to observe structure and morphology of cells in a tissue.




## EX2: Finding Cells (Click for example)






EX3: Finding  
Tissue (Click  
for example)





## EX4: Tissue Registration (Click for example)



## Resources

<https://github.com/AcejrAZ/TuPLE/tree/master/ScipyExamples>

<https://github.com/scikit-image/skimage-tutorials>