
Tucson Python Meetup

Graphical User Interface (GUI) Intro/Demo Franklin Ventura November 15, 2016

Franklin Ventura


- BS Biomedical Engineering UofA
- MS Software/Systems Engr ASU
- Ventana Medical Systems 5 yrs
 - Biomedical/Software Engineer
 - Technology and Applied Research
 - Digital Pathology
 - Cancer Diagnostics
- Applied Imaging Systems Team
 - Karl Garsha, Taras Golota








Advanced Research Needs

Stem
Flexible System










```
[acq setting(func name='set plane', func args={'plane': (0.0, 0.0), 'index': 0}),
                                                                                                                           Channels-1:
acg setting(func name='set slice', func args={'method': 'ABSOLUTE', 'z pos': 0.0}),
                                                                                                                              2-AQUA AQUA:
acq setting(func name='set channel', func args={'index': 0, 'channel': 'channel2'}),
acq setting(func name='set exposure', func args={'exposure': 20}),
                                                                                            No one
                                                                                                                              6-CY5 CY5:
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 0, 'z pos': -0.5}),
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                              1-DAPI DAPI:
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 1, 'z pos': 0.0}),
                                                                                                                              4-GOLD GOLD:
acq setting(func name='set image', func args={'binning': 2}),
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 2, 'z pos': 0.5}),
                                                                                            wants to
                                                                                                                              3-GREEN GREEN:
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                              5-RED RED:
acq setting(func name='set channel', func args={'index': 1, 'channel': 'channel4'}),
acq setting(func name='set exposure', func args={'exposure': 40}),
                                                                                                                           Image-1:
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 0, 'z pos': -0.5}),
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                              img type: OME
                                                                                            have to
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 1, 'z pos': 0.0}),
                                                                                                                              method: Reflected
acq setting(func name='set image', func args={'binning': 2}),
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 2, 'z pos': 0.5}),
                                                                                                                           Objectives-1:
acq setting(func name='set image', func args={'binning': 2}),
acq setting(func name='set channel', func args={'index': 2, 'channel': 'channel6'}),
                                                                                                                              '20.0':
                                                                                            work
acq setting(func name='set exposure', func args={'exposure': 60}),
                                                                                                                           Planes-1:
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 0, 'z pos': -0.5})
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                              - - 1056.709
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 1, 'z pos': 0.0}),
                                                                                                                                 - 1714.655
                                                                                            with this
acg setting(func name='set image', func args={'binning': 2}),
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 2, 'z pos': 0.5}),
                                                                                                                                 -28.0
acq setting(func name='set image', func args={'binning': 2}),
acq setting(func name='set channel', func args={'index': 3, 'channel': 'channell'}),
                                                                                                                              - Objectives-1
acq setting(func name='set exposure', func args={'exposure': 10}),
                                                                                                                              - FOV 001
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 0, 'z pos': -0.5}),
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                           Save Folder: C:\TEST NEW
acq_setting(func_name='set_slice', func_args={'method': 'RELATIVE', 'index': 1, 'z_pos': 0.0}),
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                           Slices-1:
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 2, 'z pos': 0.5}),
                                                                                                                              delta: 0.5
acq setting(func name='set image', func args={'binning': 2}),
acq setting(func name='set channel', func args={'index': 4, 'channel': 'channel3'}),
                                                                                                                              direction: UP
acq setting(func name='set exposure', func args={'exposure': 30}),
                                                                                            Or even
                                                                                                                              end: 0.0
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 0, 'z pos': -0.5})
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                              method: RELATIVE
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 1, 'z pos': 0.0}),
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                              slices num: 7
                                                                                           this
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 2, 'z pos': 0.5}),
                                                                                                                              start: 0.0
acq setting(func name='set image', func args={'binning': 2}),
acq setting(func name='set channel', func args={'index': 5, 'channel': 'channel5'}),
                                                                                                                           dim order:
acq setting(func name='set exposure', func args={'exposure': 50}),
                                                                                                                              Planes-1
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 0, 'z pos': -0.5}),
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                              Slices-1
acq setting (func name='set slice', func args={'method': 'RELATIVE', 'index': 1, 'z pos': 0.0}),
acq setting(func name='set image', func args={'binning': 2}),
                                                                                                                              Channels-1
acq setting(func name='set slice', func args={'method': 'RELATIVE', 'index': 2, 'z pos': 0.5}),
                                                                                                                              Image-1
acg setting(func name='set image', func args={'binning': 2})]
```

We are okay Non-programmers with this want this

```
worker = Thread(name='ImagingQueue',
                target=self.img_engine.image_creator, args=())
 worker.setDaemon(False)
worker.start()
self.system.camera.live stop()
cont, skip = True, False
for i, acq in enumerate(acq settings):
    logger.debug("Acquisition: %s Args: %s", acq.func_name,
                 acq.func args)
    getattr(self, acq.func name)(**acq.func args)
   if progress is not None:
        (cont, skip) = progress.update(i)
        if not cont or skip or cont is None:
            self.img engine.stop = True
            return False
 worker.join()
```

return self.img engine.fpath, self.img engine.image paths




Many Options for GUI development - https://wiki.python.org/moin/GuiProgramming

- Guidata
- Kivy
- O PyForms
- PyQt
- O Pywebview
- TkInter
- Traits/TraitsUI


- VTK
- WxPython
- PyJamas
- PyMUI
- PyWebkitGtk
- O PySide
- O Pyglet


What I use/used


















wxPython number wx.Choice(self, ['1', '3', '5', ...]) interval = wx.Choice(self, ['0.25', '0.5', '0.75', ...]) number text = wx.StaticText(self, "Number of Slices") interval text = wx.StaticText(self, "Interval (micron)") zrange text = wx.StaticText(self, u"Z-Range(micron)") zrange = wx.StaticText(self, "0.0")

app = wx.App(False)dialog = slices example(None) dialog.Show(True) app.MainLoop()




sizer = wx.GridBagSizer(3, 3) sizer.Add(self.number text, pos=(0,0)) sizer.Add(self.interval text, pos=(0,1)) sizer.Add(self.zrange text, pos=(0,2)) sizer.Add(self.number, pos=(1,0))

sizer.Add(self.interval, pos=(1,1)) sizer.Add(self.zrange, pos=(1,2))

TraitsUI/PyQT

```
number = Enum(1, 3, 5, ...)
interval = Enum(0.25, 0.5, ...)
zrange = Float
slice_image = Image(paths.slice)
```

dialog = slices_example()
dialog.configure_traits()



```
View(HGroup(
Item('slice image', show label=False),
VGroup(
  Item(label='Number of Slices'),
  Item('number', show label=False,
       width=0.25)),
VGroup(
  Item(label='Interval (micron)'),
  Item('interval', show label=False,
       width=0.25)),
VGroup(
  Item(label='Z-Range (micron)'),
  Item('zrange', style='readonly',
       show label=False,
       width=0.25))
```


View Connection to Model/Controller

```
wxPython
```

```
self.number.Bind(wx.EVT CHOICE, self. zrange update)
      self.interval.Bind(wx.EVT_CHOICE, self. zrange_update)
      def zrange update(self, event):
         num = int(self.number.GetStringSelection())
         inter = float(self.interval.GetStringSelection())
         zrange = ( num - 1 ) * inter
         self.zrange.SetLabel(str(zrange))
TraitsUI/PyQT
      @on trait change('number, interval')
      def _zrange_update(self):
         self.zrange = ( self.number - 1 ) * self.interval
```

Password Control

global_access = Bool(False)



Login

Item(name='calibrate_button', label='Calibrate Light',

enabled_when='access', show_label=False)

Control Calibrate This light is calibrated for inter-system reproducibility and to enable quantitative measurement of advanced, multiplexed, fluorescence assays Calibrate Light Test Stability Stability Time ("HH:MM:SS"): 0:05:00 Calibration Date: Mon, 08 Aug 2016 10:27:07

Chaco

Interactive data visualization and exploration

- Flexible drawing and layout
- Modular and extensible architecture
- Data model for ease of extension and embedding

self.image = ArrayPlotData()
self.image.set_data(self.image_data)
plot = Plot(self.image, 'top left')
plot.tools.append(PanTool(plot))
zoom = CustomZoomTool(plot, 'box')

self.image_tool =
ImageInspectorTool(plot)

im = self.acquire()
viewer.image_data = im.astype(int)

cx_Freeze - Create Executable

- All dependencies packaged together
- No Python installed

```
build_info = Executable(
    script="gui/launchpad.py",
    base='Win32GUI',
    targetName="Ventana Acquisition Engine",
    icon=r"gui\images\icon.ico")
```

```
setup(
   version=VERSION,
   name="Ventana Acquisition Engine",
   executables=[build_info])
```

Inno Setup Compiler

- Installer File
- Even easier to distribute





Thank you

- Karl Garsha
- Dustin Harshman
- Taras Golota





GUI Workshop

Tuesday November 29, 2016 6 pm Keating Room 107

